

INSTRUCTIONS FOR COLLECTING URINE FOR KRYPTOPYRROLES

Note:

Random void is required and must be collected at Wakefield Collection Centre

- Label the containers (including full name, date of birth, time and date of collection).
- Void into the foam cup provided. (Approx ½ cup is sufficient.)
- Pour into both urine containers. No more than ⅔ full (a minimum of 10 mls (1 cm) is required in both jars).
- Shake the urine containers to mix the urine with the additive.
- Wrap both containers in foil to protect from light.
- Place the specimen containers in the specimen bag provided and seal (removing air).
- Give the specimens to the collector.
- Collector to place bag into another specimen bag with ice.

Our Client Services staff will be happy to assist you with any enquiries on:

8366 2000