

INSTRUCTIONS FOR COLLECTING URINE KRYPTOPYRROLES USING A PAEDIATRIC BAG

NOTE: 2nd morning void is preferred but a subsequent specimen is acceptable.

- See enclosed instructions on how to apply the paediatric urine bag.
- Once enough urine is present, peel bag off gently and use the syringe provided to draw out the urine.
- Gently squirt the urine into the container provided (no more than half full).
- Invert the urine container a few times to mix the urine with the additive.
- Label the container (including TIME and DATE OF COLLECTION) and wrap in foil (shiny side out) to protect from light.
- Place specimen container in the specimen bag provided and seal.
- Store in the fridge until able to transport to Kent Town Collection Centre.
- Please deliver as soon as possible.
- To transport place on ice. Put a small amount of ice in the 2nd specimen bag. Then place the specimen inside this and seal.